
# **3M ESPE**

## **Elipar™ FreeLight™**

**Cordless Polymerization Unit**

**Polymériseur sans fil**

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### Safety

#### PLEASE NOTE!


Prior to installation and start-up of the unit, carefully read the instructions provided herein!

As with all technical devices, the proper function and safe operation of this unit depend on the user's compliance with the standard safety procedures as well as the specific safety recommendations presented in these Operating Instructions.

1. The unit must be used in strict accordance with the following Operating Instructions. The manufacturer accepts no liability for any damage resulting from the use of this unit for any other purpose.
2. Prior to start-up of the unit make sure that the operating voltage stated on the rating plate is compatible with the available mains voltage. Operation of the unit at a different voltage may damage the unit.
3. Irradiation of the eyes bears an inherent risk of health damage. Hence, the light must not be directed toward the eyes. Exposure must be restricted to the area of the oral cavity in which clinical treatment is intended.
4. **CAUTION!** The Elipar FreeLight unit generates high intensity light. The emitted light should be positioned directly above the material to be cured – intensive light exposure of soft tissues (gingiva, oral mucosa, and skin) should be avoided as such exposure may cause damage or irritation. If applicable, cover such areas. If exposure of soft tissues cannot be avoided, adjust the polymerization process to the light level, e.g. by shortening the polymerization times or increasing the distance between the light guide exit and the material to be cured.
5. Condensation resulting from the unit being transferred from a cold to a warm environment may be a potential risk. Hence, the unit should be turned on only after it is completely equilibrated to ambient temperature.
6. In order to avoid electric shock do not introduce any objects into the unit with the exception of replacement parts handled in accordance with the Operating Instructions.
7. Use original 3M ESPE parts exclusively to replace defective parts in accordance with these Operating Instructions. The product guarantee does not cover any damage resulting from the use of third-party replacement parts.

Table of Contents	Page
Safety	1
Product Description	3
Fields of Application	3
Technical Data	3
Charger	3
Handpiece	3
Charger and Handpiece	3
Installation of the Unit	4
Factory Settings	4
Initial Steps	4
Charger	4
Light Guide/Handpiece	4
LED Display of the Charger	4
Table Holder for the Handpiece	4
Operation	5
Selection of Exposure Time	5
Activating and Deactivating the Light	5
Positioning the Light Guide	5
Removing and Inserting the Light Guide from/into the Handpiece	5
Measurement of Light Intensity	5
Low Battery Charge Display	6
Power-Down Mode	6
Acoustical Signals	6
Operating Errors	6
Maintenance and Care	7
Replacement of the Battery	7
Handpiece/Battery Care	7
Cleaning the Light Guide	7
Cleaning the Charger and Handpiece	7
Storage of the Handpiece during Extended Periods of Non-Use	8
Disposal	8
Customer Information	8

8. Should you have any reason to suspect the safety of the unit to be compromised, the unit must be taken out of operation and labeled accordingly to prevent third parties from inadvertently using a possibly defective unit. Safety may be compromised, e.g. if the unit malfunctions or is noticeably damaged.
9. Solvents, flammable liquids, and powerful heat sources must not be kept in the vicinity of the unit as exposure to one of these factors may damage the plastic housing of the unit.
10. Do not allow cleaning agents to enter the unit as this may short-circuit the unit or cause potentially dangerous malfunction.
11. Only factory-authorized technicians are authorized to open the housing of the unit and carry out any repairs.
12. Elipar FreeLight must not be used in patients, or by users, with heart pacemaker implants who have been advised to be cautious regards their exposure to small electrical devices.
13. Do not use Elipar FreeLight in patients with a history of photobiological reactions (including individuals with urticaria solaris or erythropoietic protoporphyria) or who are currently on photosensitizing medication (including 8-methoxypsoralen or dimethylchlorotetracycline).
14. Individuals with a history of cataract surgery may be particularly sensitive to the exposure to light and should be discouraged from Elipar FreeLight treatment unless adequate safety measures, such as the use of protective goggles to remove blue light, are undertaken.
15. Individuals with a history of retinal disease should seek advice from their ophthalmologist prior to operating the unit. In operating the Elipar FreeLight unit, this group of individuals must take extreme care and comply with any and all safety precautions (including the use of suitable light-filtering safety goggles).
16. Prior to each use of the device ensure that the emitted light intensity is sufficient to safely guarantee polymerization. Use the built-in light testing area for this purpose.

## Product Description

Elipar FreeLight, manufactured by 3M ESPE, is a light source for intraoral polymerization of dental materials. The unit consists of a charger and a cordless handpiece powered by a rechargeable battery. The unit is designed for use on a table and cannot be wall-mounted.

The unit uses Light Emitting Diodes (LEDs) as the light sources. In contrast to halogen light units, the unit emits light mainly in the wavelength range of 440 to 490 nm, i. e. the relevant range for camphor quinone-containing products.

The optimal match of the wavelength range to the intended purpose ensures that the polymerization performance is similar to that of halogen light-emitting units, albeit at lower light intensity. This is associated with a desirable by-effect in that the operator is exposed to less glare.

Exposure time options:

- 10, 20, 30 or 40 sec

The charger is equipped with an integral light intensity testing area.

The device is shipped with a turbo light guide with an 8 mm diameter light exit. It is not permissible to use the light guides of other units.

The 'maxi fiber rod' with a diameter of 13 mm for larger areas, e. g. for fissure sealing, and the 'proxi fiber rod' with a point-shaped light exit, e. g. for use in interproximal areas, are available as accessories.

The handpiece is equipped with a "power-down" function to minimize the unit's energy consumption. The handpiece switches to "power-down" mode once it is placed in the charger or if left unused for approx. 10 min outside the charger.

In stand-by mode, the charger consumes maximally 0.75 W. This corresponds to the power consumption recommended by the European Commission pursuant to the "Code of Conduct" for stand-by operation (from 2003).

## Fields of Application

- Polymerization of light-curing dental materials with photo initiator for the wavelength range 440–490 nm
  - Though the majority of light-curing dental materials are responsive in this range of wavelengths, you may wish to contact the manufacturer of the material in question.

## Technical Data

### Charger

Operating voltage: 100 V, 120 V, 230 V 50/60 Hz  
(see rating plate for factory-set voltage)

Power input: max. 10 VA  
max. 0.75 W in stand-by mode

Dimensions: depth: 210 mm  
width: 95 mm  
height: 60 mm

Weight: 555 g

### Handpiece

Power supply: Nickel-metal hydride battery, 4.8 V

Wavelength range: 440–490 nm

Light intensity: approx. 400 mW/cm<sup>2</sup>

Duration of continuous use: min. 7 min (before activation of temperature control)

Total exposure time with new, fully charged battery: typically 45 min

Dimensions: diameter: 30 mm  
length: 285 mm

Weight: 220 g

### Charger and Handpiece

Time to charge empty battery: approx. 2 h

Operating temperature: 16°C ... 40°C/59°F ... 104°F

Relative humidity: max. 80% at 37°C/99°F  
max. 50% at 40°C/104°F

Total height with handpiece inserted in the charger: 190 mm

Classification: Protection class II,  Type B

*Subject to technical modification without prior notice.*

# Installation of the Unit

## Factory Settings

The factory settings of the unit are as follows:

- 40 sec exposure time

## Initial Steps


### Charger

- ▶ Please ensure first that the voltage stated on the rating plate corresponds to the existing main supply voltage. The rating plate is attached to the bottom of the unit.
- ▶ Place the charger on a level surface.
  - The vent slots on the bottom of the unit must not be obstructed to protect the device from overheating.
- ▶ Connect the power cable of the charger to the power supply.
  - The green LED on the left side of the device is illuminated. This shows that the unit is ready for operation; please refer to the section, “LED display of the charger”.

### Light Guide/Handpiece

Never insert the handpiece in the charger unless the battery is inserted in the handpiece first!


- ▶ Autoclave the light guide prior to first use.
- ▶ Then insert the light guide in the handpiece until it snaps into place.
- ▶ Insert the battery:
  - ▶ Rotate the lid at the lower end of the handpiece counterclockwise until it hits the stop, and then remove the lid.

- ▶ **Insert** the rechargeable battery in a fluent motion as directed by the arrow **until the battery can be heard to snap into place**. If it has not snapped into place properly, the rechargeable battery will cause malfunction of the unit.

- ▶ Replace the lid and rotate clockwise until it locks into place.
- ▶ In the event of malfunction, remove the rechargeable battery from the device and re-insert as described above.
- ▶ Prior to the first use, place the handpiece in the charger to fully load the new rechargeable battery.
  - The red LED of the charger illuminates after approx. 2 sec; please see also section, “LED Display of the Charger”. The battery is fully charged once the red LED is turned off.
  - Note: new batteries attain full capacity only after several charging/discharging cycles. Therefore, the battery charge may initially be sufficient for a smaller number of exposures only.

### LED Display of the Charger


LED green		LED red	
On	Off	no handpiece inserted in the charger	→ charger ready for operation
On	Off	handpiece inserted in the charger	→ charging has been completed
On	On	handpiece inserted in the charger	→ rechargeable battery is being charged
On	Flashing	handpiece inserted in the charger	→ rechargeable battery is defective

### Table Holder for the Handpiece

While performing a procedure in a patient the handpiece can be placed in a table holder.

## Operation


### Selection of Exposure Time

Exposure times of 10, 20, 30, and 40 sec are available.

- The adequate exposure time of the product at hand is stated in the product's Instructions for Use.
- Select the exposure time by pressing the button in the dark-blue field.
  - The selected exposure time is indicated by the 4 green LEDs.
  - Each time the button is pressed, the setting advances to the next (higher) value; from 40 sec the exposure time setting returns to 10 sec. Keeping the button depressed the display scrolls through the available settings.
  - While exposure is ongoing, the button for selection of the exposure time is deactivated.

### Activating and Deactivating the Light

- Activate the light by briefly pressing the Start button.
  - Initially, the LEDs show the preset exposure time: 4 illuminated LEDs represent 40 sec of exposure time. After every 10 sec of exposure, one of the LEDs is turned off, i.e. 3 LEDs correspond to 30 sec of exposure time remaining, 2 LEDs to 20 sec of exposure time remaining, etc.
- If it is desired to turn off the light before the preset exposure time has completely elapsed, this can be done by again pressing the Start button.

### Positioning the Light Guide

- Rotate the light guide into the optimal position for polymerization.
- To make full use of the light intensity provided, place the light guide as close to the filling as possible. Avoid directly contacting the filling material!
  - Keep the light guide clean at all times to obtain full light intensity.
  - Damaged light guides substantially reduce the light power and must be replaced immediately, sharp edges may cause serious injury!

### Removing and Inserting the Light Guide from/into the Handpiece

- To remove the light guide from the handpiece strongly pull towards the front of the device.
- Push the light guide into the handpiece until it snaps into position.

### Measurement of Light Intensity


The light intensity can be reliably determined only with the charger of the Elipar FreeLight unit! The light testing area is situated in the area of the charger accepting the handpiece. Because of differences in the light sources used and arrangement of the components, determination of the light intensity with the light testing components of other units yields erroneous results.

Caution!

Measure the light intensity with the turbo light guide only!

- Clean the test area with a wet cloth, if required, taking care not to **bend or damage the charging contacts!**
- Without applying pressure, place the light guide on the test area such that the light-emitting opening of the rod is level with the test area.
- Activate the lamp by pressing the Start button.
  - All 5 blue LEDs come ON for approx. 1 sec.
  - After this period of time, the number of illuminated LEDs is indicative of the measured light intensity: 5 LEDs = 100%, 4 LEDs = 80%, 3 LEDs = 60%, 2 LEDs = 40%, 1 LED = 20%.

- ▶ If the light intensity is below 80% – fewer than 4 LEDs are ON – check the light guide for contamination or defects.
- ▶ Proceed as follows:
  - (1.) Clean the light guide, if contaminated; please refer to the section on “Care” or
  - (2.) replace the light guide, if defective, or
  - (3.) if (1.) and (2.) fail to improve the situation, please contact the 3M ESPE Customer Service or your retailer.

### Low Battery Charge Display

After frequent use the battery charge of the device may drop to below 10%, so that only a few more exposures can be performed without re-charging. Low battery charge is shown by several means:

- A short audible alarm is emitted 5x at the completion of an exposure, and repeated thereafter whenever any button is touched.
  - The 10-sec-LED of the handpiece begins to flash.
- ▶ The handpiece should be placed in the charger as soon as possible to re-charge the battery.

### Power-Down Mode

Once the handpiece is placed in the charger, all internal functions and LEDs are automatically turned off as the handpiece switches to power-down mode. This reduces the power consumption of the rechargeable battery to a minimal level. Outside the charger, the handpiece is also switched to power-down mode if left unused for approx. 10 min.

- ▶ To terminate the power-down mode press one of the two buttons.
  - The power-down termination signal (two short audible signals) is emitted indicating that the

handpiece is ready for operation: the handpiece displays show the latest selected exposure mode and time settings.

### Acoustical Signals

An acoustical signal is emitted

- every time a button is pressed,
- every time the light is turned ON or OFF,
- after 10 sec of exposure time have elapsed (1 alarm signal), 2x after 20 sec, 3x after 30 sec.

Two acoustical signals are emitted

- every time the power-down mode is terminated by pressing any of the buttons.

A 2 sec-error signal is emitted, if

- the handpiece over-heats,
- the rechargeable battery is discharged to an extent that reliable polymerization is no longer ensured.

### Operating Errors

Error	Cause ▶ Solution
The 10 sec-LED within the dark blue frame flashes and a short signal is emitted 5x whenever a button is touched or the light is turned off.	The residual battery charge has dropped below 10%.  ▶ Place the handpiece in the charger and re-charge the battery.
The ongoing exposure is interrupted (light off signal is emitted) followed by a 2 sec-error signal; the device switches to power-down mode and resists further activation.	The battery lacks sufficient charge.  ▶ Place the handpiece in the charger and re-charge the battery.
A 2 sec-error signal is emitted upon pressing of the Start button.	The handpiece has become overheated in the course of use. The handpiece can be used again once it has cooled down.  ▶ Allow the handpiece to cool down. <ul style="list-style-type: none"> <li>– The handpiece can be used again, once pressing the Start button successfully activates the light.</li> </ul>

Error	Cause ► Solution
The red LED of the charger flashes.	The battery is defective. ► Replace the battery.
The green LED of the charger fails to light up, even though the power cable is connected to a power outlet.	Power outlet carries no voltage. ► Use a different power outlet. Charger is defective. ► Have the charger repaired.


## Maintenance and Care

### Replacement of the Battery


Never place the handpiece in the charger without the battery inserted in the handpiece!

Use 3M ESPE batteries only! The use of other manufacturer batteries or non-rechargeable batteries/primary batteries is a potential hazard and may damage the device.

- Rotate the lid at the lower end of the handpiece counterclockwise until it hits the stop, and then remove the lid.
- Push a small spatula or similar aid between the cylinders of the battery until the snap-lock releases so that the battery protrudes slightly from the handpiece body.

- Remove the battery from the handpiece.

- Push the (replacement) battery into the handpiece proceeding along the direction shown by the arrow until the battery locks into place.
- Replace the lid and close by rotating clockwise.
- Insert the handpiece into the unit to completely charge the new battery for its first use.
  - The red LED of the charger illuminates after approx. 2 sec; please see also “LED Display of the Charger”. The battery charging process is complete once the red LED is turned off.
  - Note: the new battery attains full capacity only after several charging/discharging cycles. Therefore, the battery charge may initially be sufficient for a smaller number of exposures only.

### Handpiece/Battery Care

- The charging contacts of the charger and handpiece must not be contacted by metallic or greasy parts.
- Do not use any other chargers as this may damage the battery!
- The battery must not be immersed in water or thrown into open fires!

### Cleaning the Light Guide

The light guide can be autoclaved.

#### Do not sterilize by chemical means or hot-air!

- The light guide should be regularly wiped clean with a soft cloth. Especially before and after steam sterilization the spots of dried liquid should be wiped off the ends of the device.
- Adhering polymerized composite should be removed with alcohol, a plastic spatula may help in removing the material.
  - Do not use any sharp or pointed tools to protect the surface of the device from scratching.

### Cleaning the Charger and Handpiece

- For disinfection of all components of the device, spray the disinfection agent onto a piece of cloth and proceed to disinfect the device.
  - Disinfection agents must not enter the device!
  - Pursept-A (from Merz) and FD 322 (from Dürre) are recommended disinfection agents.
- Clean the charger and the handpiece with a soft cloth and a mild cleaning agent, if required.
  - It is not permissible to use solvents or abrasive cleaning agents as these may damage plastic parts of the device!
  - Cleaning agents must not enter the devices!

## **Storage of the Handpiece during Extended Periods of Non-Use**

- ▶ If the handpiece is not to be used for an extended period of time – e.g. during vacation – fully charge the battery prior to departure or keep the handpiece inserted in the **operational** charger.
  - Otherwise, nearly discharged batteries may become excessively discharged due to the low degree of power consumption ongoing even in Power Down mode. This may damage the battery.
- ▶ Fully or nearly discharged batteries must be recharged as soon as possible.

## **Disposal**

As a means of protection of the environment, your new device is equipped with a Nickel-metal hydride battery. This kind of battery is free from toxic heavy metal ions.

- ▶ Dispose of defective batteries and units in accordance with local legal regulations!

## **Limited Warranty – U.S. Only**

3M ESPE guarantees your Elipar FreeLight against defects in material and workmanship for two years from date of purchase. During the two-year warranty period, 3M ESPE will repair defects or replace a defective unit free if the unit was new when it was purchased, and it is promptly returned to the following address:

### **3M Health Care Service Center**

Suite 200, Bldg. 502-1W-01  
Oakdale, MN 55128  
1-800-292-6298

### **In Canada:**

### **3M Canada Service Support Centre**

80 Enterprise Drive South  
London, Ontario N6N 1C2

Defects caused by misuse, neglect, accident, or abuse are not covered by this warranty.

This warranty does not cover lamps, fiber optic tips, and rechargeable batteries.

THIS WARRANTY IS EXPRESSLY IN LIEU OF AND EXCLUDES WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. 3M ESPE SHALL HAVE NO LIABILITY OF ANY KIND FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES OR EXPENSES THAT ARISE DIRECTLY OR INDIRECTLY

FROM THE USE OF THIS PRODUCT, OTHER THAN TO REPLACE OR REPAIR THE DEVICE. THIS LIMITED WARRANTY SHALL BE THE EXCLUSIVE REMEDY AVAILABLE TO ANY PERSON. THE FOREGOING LIMITATIONS OF REMEDIES AND LIABILITIES MAY NOT BE CHANGED EXCEPT BY WRITTEN AGREEMENT SIGNED BY AN AUTHORIZED OFFICER OF 3M ESPE.

## **Customer Information**

For questions or comments in U.S.A. or Canada please call toll-free 1-800-634-2249.

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

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3M, ESPE, Elipar and FreeLight are trademarks of 3M or 3M ESPE AG.
© 3M 2001. All rights reserved.
**3M ESPE AG**